



- Quick results: 10 s for lactate
- Microsampling sample size: 0.8 μ L
- Know your limit, be your best



Blood Lactate Meter

Eaglenos Sciences, Inc.

➤ Product features

| Easy operations & fast results



Microsampling
0.8 μ L



Quick results:
10 s for lactate



Automatic
siphon
technology



Easy operations




Easy strip ejection
no cross contamination
technology



Power-on upon
strip insertion

| Gold electrode technology



FAD-GDH
technology to
avoid blood oxygen
interference



Enzymatic biosensor
with high selectivity



Enzyme chemistry
design to prevent
ascorbic acid
interference



Sophisticated
temperature com-
pensation to ensure
accuracy and
consistency



Gold electrode,
excellent sensitiv-
ity, conductivity
and stability




Auto moisture-
detection for
damaged strip


I Accurate measurement

▪ Clinical test results of Eaglenos blood lactate test strips

Linear regression graph between Eaglenos results in fingertip whole blood and reference product results in venous whole blood



Linear regression graph between Eaglenos results in venous whole blood and reference product results in venous whole blood



Results: Lactate concentration <4.5 mmol/L, absolute difference: ± 1.2 mmol/L
 Lactate concentration ≥ 4.5 mmol/L, relative difference: $\pm 20\%$
 More than 95% of test results meet the above requirements.

I User-friendly design



Smart Identification

Automatic identification of test strips without manual adjustment



Mass Memory Storage

Secure storage up to 800 sets



Bluetooth Transmission

No manual input required, easy access to data



➤ Significance of lactate testing

Blood lactate is the final product of anaerobic glycolysis. During physical activity, as exercise intensity increases the rate of glycolysis accelerates, leading to a gradual increase in blood lactic acid production. The variation in blood lactate represents a balance between the rate of lactate production in tissues such as skeletal muscles, the rate of muscle lactate release into blood, and the rate of lactate elimination from blood. High blood lactate concentration coupled with poor performance often indicates suboptimal condition or lower aerobic metabolism capacity. Therefore, lactate testing can aid coaches in optimizing athletes' training plans, regulating athletes' exercise intensities, and assessing athletes' training levels. Hence, lactate has the reputation of being the "training guide".

➤ Parameters

Subject	Specification
Technology	Electrochemical biosensor technology
Measurement range	Blood lactate: 0.5 - 28.0 mmol/L
Sample type	Fresh fingertip capillary blood; venous whole blood
Sample volume	0.8 µL
Test time	10 s
Hematocrit range	30% - 60%
Memory	800 sets
Operating condition	15°C - 35°C RH: 10% - 85% (non-condensing) Atmospheric pressure: 70 - 106.6 kPa
Device warranty time	5 years
Auto strip ejection	✓

Contacts

Europe

☎ +41 782166266
 ✉ info.ch@eaglenos.com
 🌐 www.eaglenos.com/en
 📍 Bachtobelstrasse 5, 8810 Horgen
 Switzerland

North America

☎ +1 888-522-6398
 ✉ info@eaglenos.com
 🌐 www.eaglenos.com/en
 📍 3928 Varsity Dr., Ann Arbor, MI 48108
 USA

Asia

☎ +86 400-019-0069
 ✉ info@eaglenos.com
 🌐 www.eaglenos.com/en
 📍 B2-2, Treehouse Headquarters Cluster
 73 Tanmi Road, Nanjing, 210044, China

*Please refer to the instruction manual for further information.